

CERTIFICATE OF ANALYSIS

Complement Technology, Inc.
4801 Troup Hwy, Suite 701
Tyler, Texas 75703, USA

Product: **SGVB⁺⁺**
Catalog # **B110**
Prep. Date: _____

<u>Specifications</u>	<u>Limits</u>	<u>Results</u>
FILL VOLUME	250 – 275 mL	260 mL
BUFFER	Sucrose Gelatin Veronal Buffered saline plus calcium and magnesium, pH 7.3 ± 0.1	Sucrose Gelatin Veronal Buffered saline plus calcium and magnesium, pH 7.3 ± 0.1
PRESERVATIVE	Sodium Azide, 0.025%	Sodium Azide, 0.025%
FILTRATION	0.22 µm filter	0.22 µm filter
IONIC STRENGTH	6 mmho ± 1 @ 22 °C	6 mmho ± 1 @ 22 °C

STORAGE CONDITIONS

Buffer may be stored for 3 months at 4 °C.
Buffer may be stored for 4 years at -20 °C. If frozen, buffer should be warmed and stirred prior to use to assure that the gelatin returns completely into solution. Buffer should then be chilled to 4 °C before use.

FOR RESEARCH USE ONLY
NOT FOR HUMAN OR DRUG USE

Signature of Analyst

Date of Analysis