

CERTIFICATE OF ANALYSIS

Complement Technology, Inc.
4801 Troup Hwy, Suite 701
Tyler, Texas 75703, USA

Product: **VBS^o**
Catalog # **B114**
Prep. Date: _____

| <u>Specifications</u> | <u>Limits</u> | <u>Results</u> |
|-----------------------|---|---|
| FILL VOLUME | 250 – 275 mL | 260 mL |
| BUFFER | Veronal Buffered saline without calcium and magnesium, pH 7.3 ± 0.1 | Veronal Buffered saline without calcium and magnesium, pH 7.3 ± 0.1 |
| PRESERVATIVE | Sodium Azide, 0.025% | Sodium Azide, 0.025% |
| FILTRATION | 0.22 µm filter | 0.22 µm filter |
| IONIC STRENGTH | 14 mmho ± 1 @ 22 °C | 14 mmho ± 1 @ 22 °C |

STORAGE CONDITIONS

Buffer may be stored for 3 months at 4 °C.
Buffer may be stored for 4 years at -20 °C. If frozen, buffer should be warmed and stirred prior to use. Buffer should then be chilled to 4 °C before use.

FOR RESEARCH USE ONLY
NOT FOR HUMAN OR DRUG USE

Signature of Analyst

Date of Analysis