

CERTIFICATE OF ANALYSIS

Complement Technology, Inc.
4801 Troup Hwy, Suite 701
Tyler, Texas 75703, USA

Product: **GVB⁺⁺ No Azide**
Catalog # **B123**
Prep. Date: _____

<u>Specifications</u>	<u>Limits</u>	<u>Results</u>
FILL VOLUME	1000 – 1010 mL	1000 mL
BUFFER	Gelatin Veronal Buffered saline plus calcium and magnesium, pH 7.3 ± 0.1	Gelatin Veronal Buffered saline plus calcium and magnesium, pH 7.3 ± 0.1
PRESERVATIVE	None	None
FILTRATION	0.22 μ m filter	0.22 μ m filter
IONIC STRENGTH	14 mmho ± 1 @ 22 °C	14 mmho ± 1 @ 22 °C

STORAGE CONDITIONS

Buffer may be stored for 3 months at 4 °C.

Buffer may be stored for 4 years at -20 °C. If frozen, buffer should be warmed and stirred prior to use to assure that the gelatin returns completely into solution. Buffer should then be chilled to 4 °C before use.

**Thawing Frozen Buffer:

Thaw quickly at 37°C, continually stirring until all particles are dissolved, chilled to 4 °C.

FOR RESEARCH USE ONLY
NOT FOR HUMAN OR DRUG USE

Signature of Analyst

Date of Analysis

Phone: 1-903-581-8284 FAX: 1-903-581-0491 Email: contactCTI@Complementtech.com
Web Site: www.ComplementTech.com