

CERTIFICATE OF ANALYSIS

Complement Technology, Inc.
4801 Troup Hwy, Suite 701
Tyler, Texas 75703, USA

Product: **GVBE**
Catalog # **B104**
Prep. Date: _____

<u>Specifications</u>	<u>Limits</u>	<u>Results</u>
FILL VOLUME	250 – 275 mL	260 mL
BUFFER	Gelatin Veronal Buffered saline with 10 mM EDTA, pH 7.3 $\pm$ 0.1	Gelatin Veronal Buffered saline with 10 mM EDTA, pH 7.3 $\pm$ 0.1
PRESERVATIVE	Sodium Azide, 0.025%	Sodium Azide, 0.025%
FILTRATION	0.22 μ m filter	0.22 μ m filter
IONIC STRENGTH	14 mmho $\pm$ 1 @ 22 °C	14 mmho $\pm$ 1 @ 22 °C

STORAGE CONDITIONS

Buffer may be stored for 3 months at 4 °C.

Buffer may be stored for 4 years at -20 °C. If frozen, buffer should be warmed and stirred prior to use to assure that the gelatin returns completely into solution. Buffer should then be chilled to 4 °C before use.

****Thawing Frozen Buffer:**

Thaw quickly at 37°C, continually stirring until all particles are dissolved, chilled to 4 °C.

FOR RESEARCH USE ONLY
NOT FOR HUMAN OR DRUG USE

CERTIFICATE OF ANALYSIS

Complement Technology, Inc.

4801 Troup Hwy, Suite 701

Tyler, Texas 75703, USA

Product: **GVBE**

Catalog # **B104**

Prep. Date: _____

Signature of Analyst

Date of Analysis