

CERTIFICATE OF ANALYSIS

Complement Technology, Inc.
4801 Troup Hwy, Suite 701
Tyler, Texas 75703, USA

Product: **SGVB⁺⁺**
Catalog # **B111**
Prep. Date: _____

<u>Specifications</u>	<u>Limits</u>	<u>Results</u>
FILL VOLUME	1000 – 1100 mL	1050 mL
BUFFER	Sucrose Gelatin Veronal Buffered saline plus calcium and magnesium, pH 7.3 ± 0.1	Sucrose Gelatin Veronal Buffered saline plus calcium and magnesium, pH 7.3 ± 0.1
PRESERVATIVE	Sodium Azide, 0.025%	Sodium Azide, 0.025%
FILTRATION	0.45 µm filter	0.45 µm filter
IONIC STRENGTH	6 mmho ± 1 @ 22 °C	6 mmho ± 1 @ 22 °C

STORAGE CONDITIONS

Buffer may be stored for 3 months at 4 °C.

Buffer may be stored for 4 years at -20 °C. If frozen, buffer should be warmed and stirred prior to use to assure that the gelatin returns completely into solution. Buffer should then be chilled to 4 °C before use.

****Thawing Frozen Buffer:**

Thaw quickly at 37°C, continually stirring until all particles are dissolved, chilled to 4 °C.

**FOR RESEARCH USE ONLY
NOT FOR HUMAN OR DRUG USE**

Signature of Analyst

Date of Analysis

Phone: 1-903-581-8284 FAX: 1-903-581-0491 Email: contactCTI@aol.com
Web Site: www.ComplementTech.com